

## Frequently Asked Questions

### 1. What material is Norris Clear made from?

Injection-molded polycrystalline alumina (ceramic) for high translucency, precision, and consistent strength.

### 2. Is the door mechanism ceramic or metal?

The door mechanism is Ceramic but is the same trusted mechanism as the metal Norris door with reinforced, thicker door housing.

### 3. Will Norris Clear match my current Norris Metal prescriptions and workflows?

Yes. Rx and slot size are identical to Norris Metal and are fully compatible, brackets and wires work seamlessly together.

### 4. Can I use the same wires, ligation protocol, and accessories?

Yes. Use your existing Norris wire sequences and accessories; Norris Clear integrates into the same clinical protocol. However, it's important to let Niti wires work longer. The largest SS wire size is 17 x 25.

### 5. How do these compare esthetically and functionally to aligners?

Norris Clear provides a fully clear, fixed-appliance option with excellent esthetics and can often achieve difficult movements faster or more predictably than aligners, especially in complex cases.

### 6. What's included in starter/intro kits?

Intro kits include 10 bracket cases, Niti Infinity wires, Clear Bracket Opening Tool, HIP Tool (archwire seater), plus a complimentary Clear Debonding Plier.