

Atlas Mini Bracket System , MBT™ Rx - 10 Pack

MAXILLARY	REF	HK	TRQ	ANG	M/D	IN/OUT	.018 RIGHT	.018 LEFT	.022 RIGHT	.022 LEFT
Central	UR ● UL ●	-	17	4	3.1mm	0.40mm	#201-M8-11	#201-M8-21	#201-M2-11	#201-M2-21
Lateral	UR ● UL ●	-	10	8	2.8mm	0.78mm	#201-M8-12	#201-M8-22	#201-M2-12	#201-M2-22
Lateral w/ Hk	UR ● UL ●	D	10	8	2.8mm	0.78mm	#201-M8-12D	#201-M8-22D	#201-M2-12D	#201-M2-22D
Cuspid	UR ● UL ●	-	-7	8	2.8mm	0.40mm	#201-M8-13	#201-M8-23	#201-M2-13	#201-M2-23
Cuspid	UR ● UL ●	-	0	8	2.8mm	0.40mm	#201-M8-13C	#201-M8-23C	#201-M2-13C	#201-M2-23C
Cuspid w/ Hk	UR ● UL ●	D	-7	8	2.8mm	0.40mm	#201-M8-131	#201-M8-231	#201-M2-131	#201-M2-231
Cuspid w/ Hk	UR ● UL ●	D	0	8	2.8mm	0.40mm	#201-M8-131C	#201-M8-231C	#201-M2-131C	#201-M2-231C
1st Bicuspid	UR ● UL ●	-	-7	0	2.8mm	0.40mm	#201-M8-14	#201-M8-24	#201-M2-14	#201-M2-24
1st Bicuspid w/ Hk	UR ● UL ●	M	-7	0	2.8mm	0.40mm	#201-M8-141	#201-M8-241	#201-M2-141	#201-M2-241
1st Bicuspid DG Hk	UR ● UL ●	D	-7	0	2.8mm	0.40mm	#201-M8-141D	#201-M8-241D	#201-M2-141D	#201-M2-241D
2nd Bicuspid	UR ● UL ●	-	-7	0	2.8mm	0.40mm	#201-M8-15	#201-M8-25	#201-M2-15	#201-M2-25
2nd Bicuspid w/ Hk	UR ● UL ●	M	-7	0	2.8mm	0.40mm	#201-M8-151	#201-M8-251	#201-M2-151	#201-M2-251
2nd Bicuspid DG Hk	UR ● UL ●	D	-7	0	2.8mm	0.40mm	#201-M8-151D	#201-M8-251D	#201-M2-151D	#201-M2-251D
MANDIBULAR	REF	HK	TRQ	ANG	M/D	IN/OUT	.018 RIGHT	.018 LEFT	.022 RIGHT	.022 LEFT
Lower Anteriors	LR ● LL ●	-	-6	0	2.2mm	0.76mm	#201-M8-31	Universal	#201-M2-31	Universal
Lower Anteriors Wide	LR ● LL ●	-	-6	0	2.9mm	0.78mm	#201-M8-31S	Universal	201-M2-31S	Universal
Lateral w/ Hk	LR ● LL ●	D	-6	0	2.2mm	0.76mm	#201-M8-32D	#201-M8-42D	#201-M2-32D	#201-M2-42D
Cuspid	LR ● LL ●	-	-6	3	2.8mm	0.40mm	#201-M8-33	#201-M8-43	#201-M2-33	#201-M2-43
Cuspid w/ Hook	LR ● LL ●	D	-6	3	2.8mm	0.40mm	#201-M8-331	#201-M8-431	#201-M2-331	#201-M2-431
1st Bicuspid	LR ● LL ●	-	-12	2	2.8mm	0.40mm	#201-M8-34	#201-M8-44	#201-M2-34	#201-M2-44
1st Bicuspid w/ Hk	LR ● LL ●	M	-12	2	2.8mm	0.40mm	#201-M8-341	#201-M8-441	#201-M2-341	#201-M2-441
1st Bicuspid DG Hk	LR ● LL ●	D	-12	2	2.8mm	0.40mm	#201-M8-341D	#201-M8-441D	#201-M2-341D	#201-M2-441D
2nd Bicuspid	LR ● LL ●	-	-17	2	2.8mm	0.40mm	#201-M8-35	#201-M8-45	#201-M2-35	#201-M2-45
2nd Bicuspid w/ Hk	LR ● LL ●	M	-17	2	2.8mm	0.40mm	#201-M8-351	#201-M8-451	#201-M2-351	#201-M2-451
2nd Bicuspid DG Hk	LR ● LL ●	D	-17	2	2.8mm	0.40mm	#201-M8-351D	#201-M8-451D	#201-M2-351D	#201-M2-451D

Atlas Mini Bracket System , MBT™ Rx - 5x5 Set

FULL SETS	.018 SLOT	.022 SLOT
Upper/Lower 5x5 No Hks	#201-M8-NHK	#201-M2-NHK
Upper/Lower 5x5 Hks on 3's	#201-M8-3	#201-M2-3
Upper/Lower 5x5 Hks on 3,4,5's	#201-M8-345	#201-M2-345
Upper/Lower 5x5 Hks on 3,4,5's Distal Gingival Bicuspid	#201-M8-345DG	#201-M2-345DG

*DynaFlex® does not imply any endorsements from Drs. McLaughlin, Bennett and Trevisi and does not claim to be an exact duplicate of the MBT™ system. MBT™ is a registered trademark of 3M/Unitek Corporation.



www.dynaflex.com | 866-346-5665

Atlas Mini Bracket System With Offset Base, MBT™ Rx - 10 Pack

MAXILLARY	REF	HK	TRQ	ANG	M/D	IN/OUT	.018 RIGHT	.018 LEFT	.022 RIGHT	.022 LEFT
1st Bicuspid	UR ●	-	-7	0	2.8mm	0.40mm	#201-M8-140-2	#201-M8-240-2	#201-M2-140-2	#201-M2-240-2
	UL ●									
1st Bicuspid w/ Hk	UR ●	D	-7	0	2.8mm	0.40mm	#201-M8-1410-2	#201-M8-2410-2	#201-M2-1410-2	#201-M2-2410-2
	UL ●									
2nd Bicuspid	UR ●	-	-7	0	2.8mm	0.40mm	#201-M8-150-2	#201-M8-250-2	#201-M2-150-2	#201-M2-250-2
	UL ●									
2nd Bicuspid w/ Hk	UR ●	D	-7	0	2.8mm	0.40mm	#201-M8-1510-2	#201-M8-2510-2	#201-M2-1510-2	#201-M2-2510-2
	UL ●									
MANDIBULAR	REF	HK	TRQ	ANG	M/D	IN/OUT	.018 RIGHT	.018 LEFT	.022 RIGHT	.022 LEFT
1st Bicuspid	LR ●	-	-12	+2	2.8mm	0.40mm	#201-M8-340-2	#201-M8-440-2	#201-M2-340-2	#201-M2-440-2
	LL ●									
1st Bicuspid w/ Hk	LR ●	D	-12	+2	2.8mm	0.40mm	#201-M8-3410-2	#201-M8-4410-2	#201-M2-3410-2	#201-M2-4410-2
	LL ●									
2nd Bicuspid	LR ●	-	-17	+2	2.8mm	0.40mm	#201-M8-350-2	#201-M8-450-2	#201-M2-350-2	#201-M2-450-2
	LL ●									
2nd Bicuspid w/ Hk	LR ●	D	-17	+2	2.8mm	0.40mm	#201-M8-3510-2	#201-M8-4510-2	#201-M2-3510-2	#201-M2-4510-2
	LL ●									

Atlas Mini Bracket System With Oval Base, MBT™ Rx - 10 Pack

MAXILLARY	REF	HK	TRQ	ANG	M/D	IN/OUT	.018 RIGHT	.018 LEFT	.022 RIGHT	.022 LEFT
1st Bicuspid	UR ●	-	-7	0	2.8mm	0.40mm	#201-M8OV-14	#201-M8OV-24	#201-M2OV-14	#201-M2OV-24
	UL ●									
1st Bicuspid w/ Hook	UR ●	D	-7	0	2.8mm	0.40mm	#201-M8OV-141	#201-M8OV-241	#201-M2OV-141	#201-M2OV-241
	UL ●									
2nd Bicuspid	UR ●	-	-7	0	2.8mm	0.40mm	#201-M8OV-15	#201-M8OV-25	#201-M2OV-15	#201-M2OV-25
	UL ●									
2nd Bicuspid w/ Hook	UR ●	D	-7	0	2.8mm	0.40mm	#201-M8OV-151	#201-M8OV-251	#201-M2OV-151	#201-M2OV-251
	UL ●									
MANDIBULAR	REF	HK	TRQ	ANG	M/D	IN/OUT	.018 RIGHT	.018 LEFT	.022 RIGHT	.022 LEFT
1st Bicuspid	LR ●	-	-12	0	2.8mm	0.40mm	#201-M8OV-34	#201-M8OV-44	#201-M2OV-34	#201-M2OV-44
	LL ●									
1st Bicuspid w/ Hook	LR ●	D	-12	0	2.8mm	0.40mm	#201-M8OV-341	#201-M8OV-441	#201-M2OV-341	#201-M2OV-441
	LL ●									
2nd Bicuspid	LR ●	-	-17	0	2.8mm	0.40mm	#201-M8OV-35	#201-M8OV-45	#201-M2OV-35	#201-M2OV-45
	LL ●									
2nd Bicuspid w/ Hook	LR ●	D	-17	0	2.8mm	0.40mm	#201-M8OV-351	#201-M8OV-451	#201-M2OV-351	#201-M2OV-451
	LL ●									