

Atlas Mini Bright Gold Bracket System, Roth Rx - 10 Pack

MAXILLARY	HK	TRQ	ANG	M/D	IN/OUT	.018 RIGHT	.018 LEFT	.022 RIGHT	.022 LEFT
Central	-	+12	+5	3.4mm	0.61mm	#2194G-11	#2194G-21	#2195G-11	#2195G-21
Lateral	-	+8	+9	2.8mm	0.80mm	#2194G-12	#2194G-22	#2195G-12	#2195G-22
Cuspid w/ Hk	D	-2	+11	3.0mm	0.58mm	#2194G-131	#2194G-231	#2195G-131	#2195G-231
1st/2nd Bi w/ DG Hk	D	-7	0	2.8mm	0.59mm	#2194G-141D	#2194G-241D	#2195G-141	#2195G-241
MANDIBULAR	HK	TRQ	ANG	M/D	IN/OUT	.018 RIGHT	.018 LEFT	.022 RIGHT	.022 LEFT
Lower Anteriors	-	0	0	2.3mm	0.78mm	#2194G-31	Universal	#2195G-31	Universal
Cuspid w/ Hook	D	-11	+5	2.8mm	0.61mm	#2194G-331	#2194G-431	#2195G-331	#2195G-431
1st/2nd Bi w/ DG Hk	D	-17	0	2.8mm	0.64mm	#2194G-341D	#2194G-441D	#2195G-341	#2195G-441

Atlas Mini Bright Gold Bracket System, Roth Rx - 5x5 Sets

FULL SETS	.018 SLOT	.022 SLOT
Upper/Lower 5x5 Hks on 3,4,5's	#2194G-345	#2195G-345

Atlas Mini Bright Gold Bracket System , MBT™ Rx - 10 Pack

MAXILLARY	HK	TRQ	ANG	M/D	IN/OUT	.018 RIGHT	.018 LEFT	.022 RIGHT	.022 LEFT
Central	-	+17	+4	3.1mm	0.40mm	#2197G-11	#2197G-21	#2198G-11	#2198G-21
Lateral	-	+17	+4	3.1mm	0.78mm	#2197G-12	#2197G-22	#2198G-12	#2198G-22
Cuspid w/ Hook	D	-7	+8	2.8mm	0.40mm	#2197G-131	#2197G-231	#2198G-131	#2198G-231
1st/2nd Bi w/ DG Hk	D	-7	0	2.8mm	0.40mm	#2197G-141D	#2197G-241D	#2198G-141	#2198G-241
MANDIBULAR	HK	TRQ	ANG	M/D	IN/OUT	.018 RIGHT	.018 LEFT	.022 RIGHT	.022 LEFT
Lower Anteriors	-	-6	0	2.2mm	0.76mm	#2197G-31	Universal	#2198G-31	Universal
Cuspid w/ Hook	D	-6	+3	2.8mm	0.40mm	#2197G-331	#2197G-431	#2198G-331	#2198G-431
1st/2nd Bi w/ DG Hk	D	-12	+2	2.8mm	0.40mm	#2197G-341D	#2197G-441D	#2198G-341	#2198G-441

Atlas Mini Bright Gold Bracket System , MBT™ Rx - 5x5 Set

FULL SETS	.018 SLOT	.022 SLOT
Upper/Lower 5x5 Hks on 3,4,5's	#2197G-345	#2198G-345